

URINE

CHLAMYDIA / GONORRHOEA / TRICHOMONAS (First Catch)

The following information is provided as general guidance only.

Any questions should be directed to your healthcare provider.


Preparation:

The Kit contains one pipette and a yellow-cap molecular tube.

Do not urinate for at least 1 hour before collecting this sample.

Collection Procedure:

1. Collect only the first 10-30 ml of urine passed into the collection dish provided
2. Remove the screw cap from the yellow tube, always keeping the tube upright.
3. Press and hold the top of the pipette. Place the pipette into the urine and slowly release the top to draw urine into the pipette.
4. Place tip of the pipette into yellow tube and depress top to allow urine to flow into tube.
5. Add sample to the tube until it reaches a level between the indicated minimum and maximum fill lies.
6. Securely tighten the cap and gently invert the tube 3–4 times to mix well.
7. Pour any remaining urine into the toilet. Dispose of the collection dish and pipette in the rubbish
8. Label the tube with your details. **Last name, First name full and DOB.**
9. Place the yellow tube into the biohazard bag and the laboratory request form into the outside pocket of the bag.
10. Transport the specimen to the laboratory as soon as possible.



Please refer to our website for additional information

www.pathlab.co.nz

or alternatively contact
07 858 0795